


# Activación electroquímica de plataformas de oro para la detección de pesticidas mediante espectroelectroquímica Raman

David Ibáñez\*, María Begoña González-García, David Hernández-Santos, Pablo Fanjul-Bolado.

Metrohm DropSens S.L., Vivero de Ciencias de la Salud, C/Colegio Santo Domingo de Guzmán s/n, 33010 Oviedo (Spain).


\*[david.ibanez@metrohm.com](mailto:david.ibanez@metrohm.com)



# Introduction

Raman spectroscopy + Electrochemistry


## RAMAN SPECTROELECTROCHEMISTRY



- Raman spectroelectrochemistry combines the advantages of electrochemistry and Raman spectroscopy.
- Fingerprint technique: optical identification of reaction intermediates, unstable species or product structures, during an electrochemical measurement.
- **Dynamic and complete analysis** of electron transfer processes and complex redox reactions.
- Study of kinetics and mechanisms.
- **Structural properties** of materials.
- Redox chemistry of inorganic, organic and biological molecules.
- ...

# Raman Spectroelectrochemistry: SERS effect

Surface-enhanced Raman scattering (SERS) effect consists of the huge enhancement of the Raman intensity.




SERS band associated with the analyte are not observed (only two bands related to the working electrode)




SERS band associated with the analyte are detected  
SERS effect due to presence of metal nanostructures

Activation of metallic screen-printed electrodes (SPEs) by electrochemical routes leads to the generation of reproducible structures with excellent SERS properties




# Raman Spectroelectrochemistry: EC-SERS effect

EC-SERS effect enhanced or initiated by electrochemistry.




- Interaction of species with the surface can be potential-dependent



Increases analyte adsorption  
Minimizes non-specific adsorption


- In-situ* generation of fresh nanostructures in presence of the analyte

**This work** →



Fresh nanostructures:  
available for analyte interaction






## Set-up and instrument


# Spectroelectrochemical setup

## Detached set-up



# Spectroelectrochemical setup

## Fully integrated set-up






# Spectroelectrochemical instrument: SPELEC RAMAN

Fully integrated (SPELEC RAMAN) instrument for Raman spectroelectrochemistry




Main Specifications	
Laser Wavelength	785 nm
Laser optical power out	Up to 500 mW
Spectrometer Wavelength range	785 – 1010 nm
Spectrometer Resolution	$< 4 \text{ cm}^{-1}$
Raman shift	$0 - 2850 \text{ cm}^{-1}$
Potential range	$\pm 4 \text{ V}$
Maximum measurable current	$\pm 40 \text{ mA}$
Dimensions	25 x 24 x 11 cm (L x W x H)
Weight	3.95 Kg



## Results


Detection of three different pesticides is analyzed:

## Thiram:




- Dithiocarbamate fungicide
- Control of fungal diseases
- Protection of fruits and vegetables

## Imidacloprid



- Chloronicotinyl insecticide
- Control of pests (aphids, scale insects, whiteflies or leafhoppers)

## Chlorpyrifos




- Organophosphate insecticide
- Control of a variety of insect pests (mosquitoes, fire ants, cockroaches, horn flies, lice, fleas and turf)

## Thiram detection:



Initially, an EC-SERS methodology is applied. It consists of two steps in a single experiment:

1. Electrochemical activation of gold surface:

- Gold screen-printed electrode (DRP-220BT)
- HCl 0.1 M
- Cyclic voltammetry: +0.70 V to +1.40 V and back to -0.20 V. Scan rate 0.05 V/s



Electrochemical  
activation



AuNPs are deposited on the microstructures observed in the initial electrode


EC-SERS  
effect

## Thiram detection:

Initially, an EC-SERS methodology is applied. It consist of **two steps** in a single experiment:

### 2. Thiram detection:

- Integration time 2000 ms → 32 spectra during the whole experiment.




EC-SERS band (cm <sup>-1</sup> )	Assignment
349	S=CS deformation and CSS deformation
437	CH <sub>3</sub> NC deformation
548	C=S stretching
863	CH <sub>3</sub> N stretching
925	CH <sub>3</sub> N stretching and C=S stretching
1075	C=S stretching
1138	CN stretching and CH <sub>3</sub> rocking
1380	CN stretching and CH <sub>3</sub> deformation
1435	CN stretching and CH <sub>3</sub> rocking
1507	CN stretching and CH <sub>3</sub> rocking

## Thiram detection:

Initially, an EC-SERS methodology is applied. It consist of **two steps in a single experiment**:

### 2. Thiram detection:

- Integration time 2000 ms → 32 spectra during the whole experiment.



2.4 mg·L<sup>-1</sup> thiram in 0.1 M HCl solution

- Detection limit: 0.240 mg·L<sup>-1</sup> thiram
- ECM for Thiram in Water (US): 3 ug·L<sup>-1</sup>
- Maximum Residue Limits (MRLs) EU: from 0.1 to 10 mg·L<sup>-1</sup> depending on the food product



This methodology is not useful



## Thiram detection:

A new methodology is developed. It consists of three steps:

1. Preconcentration of the 60  $\mu\text{L}$  drop at 34  $^{\circ}\text{C}$  for 15 min:

→ Volume decreases to 25  $\mu\text{L}$  (HCl is concentrated until 0.24 M)



- Shorter times than 15 min do not concentrate enough the samples
- Longer times do not ensure the proper contact of the solution with WE, RE and CE.
- Higher temperature than 34  $^{\circ}\text{C}$  does not allow the exhaustive control of the volume solution due to the fast evaporation of the drop.

Preconcentration plays an **essential role** in the detection of **low concentrations**

## Thiram detection:

A new methodology is developed. It consists of **three steps**:

1. **Preconcentration** of the 60  $\mu\text{L}$  drop at **34 °C** for 15 min
2. **Electrochemical activation** of gold surface:
  - Gold screen-printed electrode (DRP-220BT)
  - HCl 0.1 M
  - Cyclic voltammetry: +0.70 V to +1.40 V and back to -0.20 V.

### 3. **Thiram detection**

Total time: 16 min approximately

- 15 min preconcentration
- 64 s  $\approx$  1 min spectroelectrochemical detection



Steps 2 and 3 in a single experiment




**Easy and quick protocol** instead of long and complex procedures of traditional methods


## Thiram detection:

Spectroelectrochemical detection:

### Preconcentration



### Without preconcentration



Cyclic voltammogram obtained in  $2.4 \text{ mg}\cdot\text{L}^{-1}$  thiram in  $0.24 \text{ M HCl}$  aqueous solution


Two new peaks are observed at **+1.02 V** and **+0.95 V** during the cathodic scan after the **preconcentration** of the sample

Voltametric profile **suggests** that these peaks are related to the **desorption processes of gold halide ions and aurous and auric complexes** generated in acid media.


## Thiram detection:

Spectroelectrochemical detection:

### Preconcentration



### Analysis of Raman spectra:



From +1.20 V to +1.40 V

- A new band is detected at 362 cm<sup>-1</sup> associated with  $\text{AuCl}_4^-$  adsorbed

At +0.80 V two bands are observed:

- At 270 cm<sup>-1</sup> for  $\text{Au-Cl}^-$  from the adsorbed  $\text{Cl}^-$
- At 323 cm<sup>-1</sup> related to  $\text{AuCl}_2^-$  adsorbed


Up to +1.20 V:

- Band at 270 cm<sup>-1</sup> disappears
- Band of  $\text{AuCl}_2^-$  increases and is slightly shifted to 333 cm<sup>-1</sup> due to the generation of different species.
- A new broad band at 570 cm<sup>-1</sup> related to **gold oxide** appears

## Thiram detection:

Spectroelectrochemical detection:

### Preconcentration




### Analysis of Raman spectra:

At lower potential than +1.00 V:


- Bands at  $333\text{ cm}^{-1}$  and  $270\text{ cm}^{-1}$  increase
- Bands at  $362\text{ cm}^{-1}$  and  $570\text{ cm}^{-1}$  disappear.

Electrodissolution of gold in acid media




## Thiram detection:

Spectroelectrochemical detection:




Comparison demonstrates that **preconcentration** of the sample **is crucial** to enhance the Raman intensity.

Raman **band at 1380 cm<sup>-1</sup>** offers the **best features** to detect lower thiram concentration due to its intensity is much higher than the other ones.

 Preconcentration methodology allows the **detection of 2.4 µg·L<sup>-1</sup> thiram** when an accurate adjustment of the baseline by polynomial fitting.



The same methodology was followed to detect **imidacloprid**:



EC-SERS band (cm <sup>-1</sup> )	Assignment
829	CCC stretching
996	CCN and CCC bending
1107	CCC stretching and CH bending
1267	NN and stretching and CH wagging
1510	CCC and CH bending
1552	ONO stretching and CCC stretching


Analysis of band at **1107 cm<sup>-1</sup>** allows the **detection of 25 μg·L<sup>-1</sup>**

ECM for Imidacloprid in Water (US) 50 ug·L<sup>-1</sup>  
MRLs (EU) from 0.05 mg·L<sup>-1</sup> depending on food



This procedure offers the **sensitivity** required.

The same methodology was followed to detect **chlorpyrifos**:




Analysis of Raman signal allows the **detection of 3.5 mg·L<sup>-1</sup>**

EC-SERS band (cm <sup>-1</sup> )	Assignment
530	P-O stretching
583	P=S and C-Cl stretching
604	P=S and C-Cl stretching
674	C-Cl stretching
749	P-O-C stretching
876	P-O-C stretching
954	Cl-ring wagging
1034	P-O-C stretching
1129	ring breathing
1264	ring vibration and CH deformation
1300	ring vibration, CH deformation and CC stretching
1406	Cl-ring vibration, CN stretching and CH deformation
1478	ring vibration, C=C stretching and CH deformation
1504	ring vibration

## Tap water analysis:


- As proof of concept, three different concentrations of thiram, 1, 3 and 20  $\mu\text{g}\cdot\text{L}^{-1}$  were prepared in 0.1 M HCl tap water solution.
- The same procedure (preconcentration and electrochemical activation) was followed:



- Analysis of band at 1380  $\text{cm}^{-1}$  shows that **1  $\mu\text{g}\cdot\text{L}^{-1}$  thiram is not detected**. On the other hand, **3 and 20  $\mu\text{g}\cdot\text{L}^{-1}$  are easily detected**.



Agreement with the previous results (2.4  $\mu\text{g}\cdot\text{L}^{-1}$  minimum value detected).



Final remarks

# Final remarks

- ✓ Raman spectroelectrochemistry based on **EC-SERS effect** offers an **easy and quick procedure** for detection of pesticides. Complicated instrumentation, tedious pretreatment protocols or time-consuming measurement procedures used in other methods are avoided.
- ✓ EC-SERS procedure proposed in this work demonstrates that the **combination of preconcentration and the electrochemical activation of gold SPEs** allows the detection of **low concentration** of pesticides with **different chemical structure**.
- ✓ Sensitivity of this methodology is demonstrated since the detection of  **$2.4 \mu\text{g}\cdot\text{L}^{-1}$  thiram and  $25 \mu\text{g}\cdot\text{L}^{-1}$  imidacloprid** is achieved.



<https://www.youtube.com/watch?v=uUkhJGFzGAI>



Thank you for you  
attention!